

ASPEN VISTA® CTO4 UPGRADE KIT

CERVICAL THORACIC ORTHOSIS UPGRADE KIT
ZERVIKAL-THORAKALE ORTHESE ERWEITERUNGSSATZ KIT DE ACTUALIZACIÓN DE ÓRTESIS CERVICAL TORÁCICA KIT DI AGGIORNAMENTO PER L'ORTESI TORACICA CERVICALE KIT DE MISE À NIVEAU DES ORTHÈSES CERVICO-THORACIQUES CERVICALE THORACALE ORTHÈSE UPGRADE KIT SADA PRO UPGRADE KRČNI HRUDNÍ ORTÉZY
OPGRADERINGSÆT TIL HALSHVIRVEL BRYSTKASSE ORTOSE ZET ANABAOMIZHĒ AYXENIKHĒ ΘΩΡΑΚΙΚΗΣ ΟΡΘΩΣΗΣ CERVIKALIT OPGRADERINGSSET FOR BRYSTORTOSE ORTEZA PIERSIOWO-SZYJNA ZESTAW ROZSZERZONY UPPGRADERINGSSATS FÖR CERVIKAL THORAKAL ORTOS КОМПЛЕКТ ДЛЯ РАШИРЕНИЯ ШЕЙНОГО КОРСЕТА KIT DE ATUALIZAÇÃO DE ÓRTESE CERVICAL TORÁCICA
颈椎畸形升级套件


ASPEN VISTA MULTIPOST BACK PANEL

ASPEN VISTA CTO4 POSTERIOR STRUT

Pat. aspenmp.com/patents
 Aspen Medical Products, LLC
 6481 Oak Canyon, Irvine, CA 92618-5202 USA
 aspenmp.com • +1-800-295-2776 • +1-949-681-0200
 ©2023 Aspen Medical Products, LLC. All Rights Reserved.
 AC0130F revised 03/2023

ENGLISH
1. STRUT APPLICATION
NOTE: If MultiPost Collar has been applied, while maintaining proper spinal motion restriction protocols, remove the posterior panel and discard. Replace the posterior panel with the one provided and attach the posterior strut (see MultiPost instructions at www.aspenmp.com).
 A. Align three connection posts with three openings in MultiPost back panel and push into position.
 B. Push down causing strut support to lock into place. Ensure that all three connectors are engaged.
 C. Slide posterior strut into slot on lumbar belt posterior panel until it snaps into place.
2. CTO 4 POSTERIOR APPLICATION
NOTE: Maintain neutral alignment, and with proper support, log roll patient onto their side.
 A. Position posterior assembly and collar back panel on patient.
 B. To adjust the length push down on button. On backside, with other hand, pry back to fully release. Adjust strut up or down.
 C. If necessary, bend metal strut to account for kyphosis.
NOTE: Confirm adjustment buttons are secure following any adjustment.
 D. Log roll patient back onto their back.
 E. Tighten both collar straps as detailed in Vista MultiPost instructions.
 F. Snap lumbar straps into vest and tighten as needed to achieve a snug, comfortable fit.
NOTE: For Vista CTO 2 application, see Vista CTO 2 instructions www.aspenmp.com.

INTENDED USE
 To provide motion restriction of the cervical spine.

INDICATIONS
 Decompressive procedures, Fracture management, Postsurgical stabilization, Whiplash, Sprains/strains, Cervicogenic headache, Cervical disc syndrome, Musculoskeletal pain, Radiculopathy.

CONTRAINDICATIONS
 Hypersensitivity or allergies to any of the materials from which the brace is made. Swelling of the lymphatic tissues caused by circulatory disorders. Patients who are not allowed to wear a brace according to medical instructions.

PRECAUTIONS
READ INSTRUCTIONS BEFORE USE. Proper training in the use of this device should take place before it is applied. These directions are guidelines only and are not offered as medical recommendations. If you suffer from a serious medical condition, we strongly suggest that you consult with a licensed health care professional before using this product. Proper fitting is required for this product to be effective. Under some circumstances, this product may be prescribed by a physician. Please see the limited warranty for further information.

CARE INSTRUCTIONS
 Hand wash only at 30° C; Wash with mild detergent; Air dry; Do not use bleach or other chemicals.

MATERIALS
 ABS, Nylon, HDPE, Polycarbonate, Delrin, Polyurethane, Cotton, Aluminium.

NOT MADE WITH NATURAL RUBBER LATEX.

DISPOSAL
 Product may be safely disposed in accordance with local laws.

REPORTING
 Please report any serious incident involving the use of this device to both the Competent Authority in your state and to the manufacturer (listed on this document).

LIMITED WARRANTY
 Aspen Medical Products, LLC, Irvine, CA 92618, warrants to the user who originally purchases this product that it is free from defects in material and workmanship. The sole obligation of Aspen Medical Products, LLC in the event of breach of warranty shall be to repair or replace the defective product or parts).

Aspen Medical Products, LLC shall have no obligation under this limited warranty in the event:
 (a) The product was not purchased from Aspen Medical Products, LLC or through its authorized channels of distribution;
 (b) The product is altered;
 (c) Any parts not supplied by Aspen Medical Products, LLC are inserted into the product; or
 (d) The product is not used in accordance with the Aspen Medical Products, LLC Instructions for Use.

THE FOREGOING IS THE SOLE AND EXCLUSIVE REMEDY FOR FAILURE IN SERVICE OF, OR DEFECT IN, THE PRODUCT. ASPEN MEDICAL PRODUCTS, LLC SHALL NOT BE LIABLE UNDER THIS OR ANY IMPLIED WARRANTY FOR ANY DIRECT, SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES. THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING THE WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR USE, AND ALL OBLIGATIONS OR LIABILITIES ON THE PART OF ASPEN MEDICAL PRODUCTS FOR DAMAGES ARISING OUT OF OR IN CONNECTION WITH THE USE OF THE PRODUCTS, WHICH ARE HEREBY DISCLAIMED AND EXCLUDED BY ASPEN MEDICAL PRODUCTS, LLC.

This warranty gives you specific legal rights and you may have other rights which vary from location to location.

A listing of authorized distributors in your area is available upon request.

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