




ASPEN EVERGREEN™ 621 SI BELT

SACRAL ORTHOSIS
SAKRALE ORTHESE
ÓRTESES SACRA
ORTESI SACRALE
ORTHÈSES SACRÉE
HEILIGBEEN ORTHESE
SAKRÁLŇÍ ORTÉZA
SAKRAL ORTOSÉ
ΟΡΘΩΣΗ ΙΕΡΟΥ ΟΣΤΟΥ
KORSRYGGORTROSE
ORTEZA KRZYŻOWA
SAKRAALINEN ORTOOSI

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ENGLISH
1. PLACEMENT
 A. Position the brace approximately 1" below the widest part of the hips (SI Joint). Center the two back panels and position symmetrically on the back.
 B. Overlap the ends of the brace and fasten by pressing down end to secure brace.
2. TIGHTEN
 To apply compression, grasp the pull tabs with the thumbs and forefingers.
3. ADJUST
 After desired compression is achieved, fasten the pull tabs to the loop material on the side panels to secure.

| SIZING GUIDE | | |
|--------------|---------------|--------------|
| Size | Circumference | |
| SMALL | 25 – 41 in | 64 – 104 cm |
| MEDIUM | 35 – 50 in | 89 – 127 cm |
| LARGE | 45 – 60 in | 114 – 152 cm |

INTENDED USE
 To provide sacral stabilization.
INDICATIONS
 Symphysis pubis dysfunction, Sacroiliac joint dysfunction, Chronic sacroiliac support, Low back sprains/strains, Sacroiliac joint fusion, Piriformis syndrome, Pelvic muscle pain, Sacroiliac pain, Low back pain, Sacroiliitis, Arthritis, Sciatica.

CONTRAINDICATIONS
 Hypersensitivity or allergies to any of the materials from which the brace is made. Swelling of the lymphatic tissue caused by circulatory disorders. Patients who are not allowed to wear a brace according to medical instructions.

PRECAUTIONS
 READ INSTRUCTIONS BEFORE USE. Proper training in the use of this device should take place before it is applied. These directions are guidelines only and are not offered as medical recommendations. If you suffer from a serious medical condition, we strongly suggest that you consult with a licensed health care professional before using this product. Proper fitting is required for this product to be effective. Under some circumstances, this product may be prescribed by a physician. Please see the limited warranty for further information.

CARE INSTRUCTIONS
 Hand wash only at 30° C; Wash with mild detergent; Air dry; Do not use bleach or other chemicals.

MATERIALS
 Nylon, Polyester, HDPE, PVC, Polyurethane, UHMWPE.

NOT MADE WITH NATURAL RUBBER LATEX.

DISPOSAL
 Product may be safely disposed in accordance with local laws.

REPORTING
 Please report any serious incident involving the use of this device to both the Competent Authority in your state and to the manufacturer (listed on this document).

LIMITED WARRANTY
 Aspen Medical Products, LLC, Irvine, CA 92618, warrants to the user who originally purchases this product that it is free from defects in material and workmanship. The sole obligation of Aspen Medical Products, LLC in the event of breach of warranty shall be to repair or replace the defective product or parts(s).

Aspen Medical Products, LLC shall have no obligation under this limited warranty in the event:
 (a) The product was not purchased from Aspen Medical Products, LLC or through its authorized channels of distribution;
 (b) The product is altered;
 (c) Any parts not supplied by Aspen Medical Products, LLC are inserted into the product; or
 (d) The product is not used in accordance with the Aspen Medical Products, LLC Instructions for Use.

THE FOREGOING IS THE SOLE AND EXCLUSIVE REMEDY FOR FAILURE IN SERVICE OF, OR DEFECT IN, THE PRODUCT. ASPEN MEDICAL PRODUCTS, LLC SHALL NOT BE LIABLE UNDER THIS OR ANY IMPLIED WARRANTY FOR ANY DIRECT, SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES. THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING THE WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR USE, AND ALL OBLIGATIONS OR LIABILITIES ON THE PART OF ASPEN MEDICAL PRODUCTS FOR DAMAGES ARISING OUT OF OR IN CONNECTION WITH THE USE OF THE PRODUCTS, WHICH ARE HEREBY DISCLAIMED AND EXCLUDED BY ASPEN MEDICAL PRODUCTS, LLC.

This warranty gives you specific legal rights and you may have other rights which vary from location to location.

A listing of authorized distributors in your area is available upon request.

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