

## ASPEN EVERGREEN™ 456 TLSO ADJUSTABLE

**THORACIC LUMBAR SACRAL ORTHOSIS  
LUMBAL-SAKRAL-THORAKALE ORTHESE  
ÓRTE SIS TORÁCICA LUMBOSACRA  
ORTESI TORÁCICA SACRO-LOMBARE  
OR THÈSE THORACIQUE LOMBO-SACRÉE  
LUMBALE SAKRALE THORACALE ORTHÈSE  
BEDERNÍ KRÍŽOVÁ HRUDNÍ ORTÉZA  
LEÑDENHVIRVELSØJLEN THORAXORTOSE  
ΟΣΦΥΓΙΚΗ ΘΩΡΑΚΙΚΗ ΟΡΘΩΣΗ ΙΕΡΟΥ ΟΣΤΕΟΥ  
LUMBAL BRYSTORTOSE  
ORTEZA PIERSIOWO-LĘDŹYWIOWO-KRZYŻOWA  
LANNERANGAN SAKRAALINEN RINTAKEHÄÄRTÖÖSI  
LUMBAL SAKRAL THORAKAL ORTOS  
ОРТЕЗА ДЛЯ ГРУДНОГО И ПОЯСНИЧНО-КРЕСТЦОВОГО ОТДЕЛОВ  
ПОВЗОЖНИКА  
ÓRTESE TORÁCICA SACRO LOMBAR  
腰骶胸廓矫形器**

Pat. aspenmp.com/patents  
Aspen Medical Products, LLC  
6481 Oak Canyon, Irvine, CA 92618-5202 USA  
aspenmp.com • +1-800-295-2776 • +1-949-681-0020  
©2021 Aspen Medical Products, LLC. All Rights Reserved.  
EG0077D revised 04/2021

**ENGLISH**  
**1. SIZE ADJUSTMENT**  
To adjust the size of the Evergreen, peel up the inside end of the side panel on each side. Adjust to achieve a proper fit and reattach.

SIZING GUIDE			
Size		Circumference	
1	28 – 37 in	71 – 94 cm	
2	36 – 47 in	91 – 119 cm	
3	46 – 54 in	117 – 137 cm	

**2. PREPARATION**  
Place the pull tabs 2 inches (5 cm) from the plastic in the back of the Evergreen brace. Pull the ends of the Evergreen brace away from each other ensuring that it is fully extended.

**3. APPLICATION**  
A. Apply the fully extended brace making sure the back panel is centered and low on the back. Apply the brace tightly, with the right side overlapping the left.  
B. When applied properly, both the upper and lower edge of the brace should be snug. If they are not, reapply the brace. If needed, angle the ends to insure a snug fit at both the upper and lower edges.  
**NOTE:** For certain patients, the ends must be steeply angled.

**4. TIGHTENING**  
A. Slip your thumbs through the holes in the pull tabs. Pull away from your body.  
B. The left pull tab controls lower compression and the right pull tab controls the upper compression. Adjust the compression as needed.  
C. When tightened, secure the pull tabs to the hook-and-loop in front of the brace.

**5. ADJUSTMENT**  
A. Adjust the upper and/or lower compression by loosening or tightening the pull tabs as needed.  
B. When patient is seated, slightly loosening the upper portion may be beneficial.

**6. STRAP ADJUSTMENT**  
A. Properly fit shoulder straps.  
B. Pull and tighten the shoulder straps.

**7. PANEL REMOVAL**  
If necessary, remove back panel after disengaging hook and loop attachment between back panel and belt.

**8. MODIFICATION OPTIONS**  
This product may be substantially modified by an individual with expertise. For more information, please visit [www.aspenmp.com/education/product-modifications](http://www.aspenmp.com/education/product-modifications).

**INTENDED USE**  
To provide motion restriction and stability of the Thoracic-lumbar-sacral spine.

**INDICATIONS**  
Postoperative support, Degenerative disc disease, Bulging or herniated disc, Fracture management, Spondylolysis, Spondylolisthesis, Spinal stenosis, Kyphosis, Spondylolysis (Osteoarthritis/Degenerative Joint Disease).

**CONTRAINDICATIONS**  
Hypersensitivity or allergies to any of the materials from which the brace is made. Swelling of the lymphatic tissue caused by circulatory disorders. Patients who are not allowed to wear a brace according to medical instructions.

**PRECAUTIONS**  
READ INSTRUCTIONS BEFORE USE. Proper training in the use of this device should take place before it is applied. These directions are guidelines only and are not offered as medical recommendations. If you suffer from a serious medical condition, we strongly suggest that you consult with a licensed health care professional before using this product. Proper fitting is required for this product to be effective. Under some circumstances, this product may be prescribed by physician. Please see the limited warranty for further information.

**CARE INSTRUCTIONS**  
Hand wash only at 30° C. Wash with mild detergent; Air dry; Do not use bleach or other chemicals.

**MATERIALS**  
Nylon, Polyester, HDPE, Polyurethane, UHMWPE.

**NOT MADE WITH NATURAL RUBBER LATEX.**

**DISPOSAL**  
Product may be safely disposed in accordance with local laws.

**REPORTING**  
Please report any serious incident involving the use of this device to both the Competent Authority in your state and to the manufacturer (listed on this document).

**LIMITED WARRANTY**  
Aspen Medical Products, LLC, Irvine, CA 92618, warrants to the user who originally purchases this product that it is free from defects in material and workmanship. The sole obligation of Aspen Medical Products, LLC in the event of breach of warranty shall be to repair or replace the defective product or parts).

Aspen Medical Products, LLC shall have no obligation under this limited warranty in the event: (a) The product was not purchased from Aspen Medical Products, LLC or through its authorized channels of distribution; (b) The product is altered; (c) Any parts not supplied by Aspen Medical Products, LLC are inserted into the product; or (d) The product is not used in accordance with the Aspen Medical Products, LLC Instructions for Use.

THE FOREGOING IS THE SOLE AND EXCLUSIVE REMEDY FOR FAILURE IN SERVICE, OR DEFECT IN, THE PRODUCT ASPEN MEDICAL PRODUCTS, LLC SHALL NOT BE LIABLE UNDER THIS OR ANY IMPLIED WARRANTY FOR ANY DIRECT, SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES. THIS WARRANTY IS VOID FOR ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING THE WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR USE, AND ALL OBLIGATIONS OR LIABILITIES ON THE PART OF ASPEN MEDICAL PRODUCTS FOR DAMAGES ARISING OUT OF OR IN CONNECTION WITH THE USE OF THE PRODUCTS, WHICH ARE HEREBY DISCLAIMED AND EXCLUDED BY ASPEN MEDICAL PRODUCTS, LLC.

This warranty gives you specific legal rights and you may have other rights which vary from location to location.

A listing of authorized distributors in your area is available upon request.

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