HORIZON[™] PRO 627

LUMBAR

DOCUMENTATION WORKSHEET: RETAIN IN PATIENT RECORD

Doctor:	Fitter:		
Patient Name:	Date:		
Patient #:	Additional Follow-Up Dates:		

Page 1 of 2


TOOLS NECESSARY: Scissors • Heat Gun • Tape Measure

FOR USE WITH PRODUCTS MANUFACTURED BY ASPEN MEDICAL PRODUCTS, LLC ONLY. THIS PRODUCT IS INTENDED FOR APPLICATION BY HEALTH CARE PRACTITIONERS AS DIRECTED BY A PHYSICIAN OR OTHER QUALIFIED MEDICAL AUTHORITY. THIS IS A PREFABRICATED ORTHOSIS. IT IS INTENDED TO BE CUSTOMIZED TO AN INDIVIDUAL PATIENT. FOLLOW THE STEPS BELOW TO CUSTOMIZE.

STEP 1 - MEASUREMENTS		STEP 2 - CUSTOMIZE BACK PANEL TO ANATOMY
1 Lower rib circumference =		A. Measure patient's lordosis then customize back panel to anatomy. B. To customize back panel, remove the panel, heat, trim, and reassemble. FLAT SIDE Patient's Lordosis
3 L1 - L5 =		Degree: ——
TIME SPENT:	2	Heat form to individual patient's anatomy and contour to create intimate fit for individual lordosis and soft tissue. Trim for individual patient's anatomy based on 3
		TIME SPENT:
STEP 3 - CUSTOMIZE SIZING	AND TIGHTENING MECHANISM	SIZING IS CRITICAL TO PROPER PERFORMANCE Use the measurements below to customize to patient's anatomy.

A. ____

B. Adjust belt to chalk mark.



C.	Adjust length of tightening mechanism. For individual patient, it
	may be necessary to adjust length of closure string. Trim and
	adjust length of strings.
	YES. AMOUNT CUT
	□ NO

TIME SPENT:	



HORIZON[™] PRO 627

LUMBAR

DOCUMENTATION WORKSHEET: RETAIN IN PATIENT RECORD

Doctor:	Fitter:		
Patient Name:	Date:		
Patient #:	Additional Follow-Up Dates:		

TOOLS NECESSARY: Scissors • Heat Gun • Tape Measure

STEP 4 - MODIFY RIGID PANELS

MODIFY ANTERIOR PANEL AS NECESSARY



TIME SPENT:				
Remove and heat mold anterior panel as necessary.				
Remove and trim to accommodate small and extra small anatomy.	9			


Page 2 of 2

STEP 5 - CUSTOMIZE BELT FIT


ANGLE ANTERIOR PANELS

Every patient has a unique individual anatomy. Determine angulation for proper fit. Circumferential contact at both upper and lower margins of brace is essential for proper brace performance and support.

- A. Bend anterior panel to conform to patient's anatomy.
- **B.** Angle anterior panels:




Configuration for best support



Inferior Angulation

Configuration for best support

TIME SPENT: _____



Superior Angulation

Configuration for best support

STEP 6 - EDUCATION

EDUCATE PATIENTS

 $Proper\ education\ is\ needed\ for\ individual\ to\ maintain\ proper\ fit\ throughout\ total\ time\ of\ wear.$

Items to educate patients on:

Independent compression mechanics

Don and doffing

	Proper	angulation	to	ensure	circumferential	contact
--	--------	------------	----	--------	-----------------	---------

Proper placement of brace

L		Proper	cleaning
---	--	--------	----------

Follow up appointments

TIME SPENT:

CLINICAL JUSTIFICATION FOR CUSTOMIZING BRACE

TOTAL TIME TO CUSTOMIZE BRACE:

