

Doctor: _____ Fitter: _____

Patient Name: _____ Date: _____

Patient #: _____ Additional Follow-Up Dates: _____

TOOLS NECESSARY: Scissors • Heat Gun • Tape Measure

FOR USE WITH PRODUCTS MANUFACTURED BY ASPEN MEDICAL PRODUCTS ONLY. THIS PRODUCT IS INTENDED FOR APPLICATION BY HEALTH CARE PRACTITIONERS AS DIRECTED BY A PHYSICIAN OR OTHER QUALIFIED MEDICAL AUTHORITY. THIS IS A PREFABRICATED ORTHOSIS. IT IS INTENDED TO BE CUSTOMIZED TO AN INDIVIDUAL PATIENT. FOLLOW THE STEPS BELOW TO CUSTOMIZE.


STEP 1 - MEASUREMENTS

1 Lower rib circumference = _____

2 Hip circumference = _____

3 L1 - L5 = _____


TIME SPENT: _____




STEP 2 - CUSTOMIZE PANELS TO ANATOMY

- A.** Measure patient's lordosis then customize panels to anatomy.
- B.** To customize panels, remove the panel, heat, trim, and reassemble.

POSTERIOR PANEL



ANTERIOR PANEL



Patient's Lordosis Degree: _____

Heat form to individual patient's anatomy and contour to create intimate fit for individual lordosis and soft tissue. Trim for individual patient's anatomy based on **3** _____

TIME SPENT: _____

STEP 3 - CUSTOMIZE SIZING AND TIGHTENING MECHANISM

SIZING IS CRITICAL TO PROPER PERFORMANCE
Use the measurements below to customize to patient's anatomy.


- A.** Use waist circumference (average of **1** and **2** _____) to determine size of Summit brace.

X-SMALL 21-27 in 46-61 cm	SMALL 26-32 in 68-81 cm	MEDIUM 31-37 in 79-94 cm
LARGE 36-42 in 91-107 cm	X-LARGE 41-47 in 104-119 cm	XX-LARGE 46-57 in 117-145 cm

- B.** Adjust length of tightening mechanism. For individual patient, it may be necessary to adjust length of closure string. Trim and adjust length of strings.

YES. AMOUNT CUT _____

NO



A. _____

TIME SPENT: _____

Doctor: _____ Fitter: _____

Patient Name: _____ Date: _____

Patient #: _____ Additional Follow-Up Dates: _____


TOOLS NECESSARY: Scissors • Heat Gun • Tape Measure

STEP 4 - CUSTOMIZE BELT FIT

ANGLE ANTERIOR PANELS


Every patient has a unique individual anatomy. Determine angulation for proper fit. Circumferential contact at both upper and lower margins of brace is essential for proper brace performance and support.

- A. Bend anterior panel to conform to patient's anatomy.
- B. Angle anterior panels:




Neutral

Configuration for best support



Inferior Angulation

Configuration for best support



Superior Angulation

Configuration for best support

TIME SPENT: _____

STEP 5 - EDUCATION

EDUCATE PATIENTS

Proper education is needed for individual to maintain proper fit throughout total time of wear.

Items to educate patients on:

- | | | |
|--|--|---|
| <input type="checkbox"/> Independent compression mechanics | <input type="checkbox"/> Proper angulation to ensure circumferential contact | <input type="checkbox"/> Proper cleaning |
| <input type="checkbox"/> Don and doffing | <input type="checkbox"/> Proper placement of brace | <input type="checkbox"/> Follow up appointments |

TIME SPENT: _____

TOTAL TIME TO CUSTOMIZE BRACE: _____