

ASPEN VISTA® CERVICAL COLLAR & ASPEN VISTA® TX COLLAR

- CERVICAL COLLAR
- ZERVIKALORTHESE
- CUELLO CERVICAL
- COLLARE CERVICALE
- COLLIER CERVICAL
- NEKBRACE
- KRČŇŮLÍMEC
- HALSHVIRVELKRAVE
- АУХЕНИКО КОЛАРО
- CERVICAL KRAGE
- SZYJNY KOLNIERZ ORTOPEDYCZNY
- KAULATUKI
- NACKKRAGE
- ШЕЙНИЙ КОРСЕТ
- COLLAR CERVICAL
- 頸托

ASPEN VISTA CERVICAL COLLAR

ASPEN VISTA TX COLLAR

MEDICAL DEVICE

SINGLE PATIENT - MULTIPLE USE

CAUTION





RESTRICTED USE

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands

Pat. aspenmp.com/patents

Aspen Medical Products, LLC
6481 Oak Canyon, Irvine, CA 92618-5202 USA
aspenmp.com • +1-800-295-2776 • +1-949-681-0200
©2021 Aspen Medical Products, LLC. All Rights Reserved.
AC0028G revised 04/2021

ENGLISH

1. PREPARATION
With the head in neutral alignment, press the back panel down onto the bed and slide it behind the patient's neck. The loop straps should be between the patient's ear and the trapezius.

2. PLACEMENT
The front panel should be at the lowest setting to begin. While pulling the sides of the front panel apart, position the bottom of the tracheal aperture at the sternal notch. Position the sides up and over the trapezius muscles. It is OK if the chin piece is not touching the chin.

3. SIZING
While holding the collar against the chest, pull out the dial to unlock. Turn the dial clockwise to raise chin piece. When the chin piece is supporting the chin, release the dial and it will self-lock.

4. TIGHTENING
A. While holding the front panel in place, attach the loop strap on each side.
B. To tighten, anchor your fingers in on the collar and peel back the loop strap on one side, tighten and reattach. Tighten the other side equally and repeat as needed to ensure a snug, symmetrical fit.
C. Check to see that the back of the chin piece is not pressing inward on the throat. If it is, lift the back ends of the side panels up and off the trapezius muscles.

TIPS

- 1. The patient's chin should be close to the front of the chin piece but the back of the chin piece should not press on the patient's throat.
- 2. All slack between the front and back panels must be removed to ensure a secure fit.
- 3. The patient often changes position, so it may be desirable to readjust the height and the back panel straps.

INTENDED USE

To provide motion restriction of the cervical spine.

INDICATIONS

Decompressive procedures, Fracture management, Postsurgical stabilization, Whiplash, Sprains/strains, Cervicogenic headache, Cervical disc syndrome, Musculoskeletal pain, Radiculopathy.

CONTRAINDICATIONS

Hypersensitivity or allergies to any of the materials from which the brace is made. Swelling of the lymphatic tissue caused by circulatory disorders. Patients who are not allowed to wear a brace according to medical instructions.

PRECAUTIONS

READ INSTRUCTIONS BEFORE USE. Proper training in the use of this device should take place before it is applied. These directions are guidelines only and are not offered as medical recommendations. If you suffer from a serious medical condition, we strongly suggest that you consult with a licensed health care professional before using this product. Proper fitting is required for this product to be effective. Under some circumstances, this product may be prescribed by a physician. Please see the limited warranty for further information.

CARE INSTRUCTIONS

Hand wash only at 30° C. Wash with mild detergent; Air dry; Do not use bleach or other chemicals.

MATERIALS

ABS, Nylon, HDPE, Polycarbonate, Delrin, Polyurethane, Cotton.

NOT MADE WITH NATURAL RUBBER LATEX.

DISPOSAL

Product may be safely disposed in accordance with local laws.

REPORTING

Please report any serious incident involving the use of this device to both the Competent Authority in your state and to the manufacturer (listed on this document).

LIMITED WARRANTY

Aspen Medical Products, LLC, Irvine, CA 92618, warrants to the user who originally purchases this product that it is free from defects in material and workmanship. The sole obligation of Aspen Medical Products, LLC in the event of breach of warranty shall be to repair or replace the defective product or parts).

Aspen Medical Products, LLC shall have no obligation under this limited warranty in the event:
(a) The product was not purchased from Aspen Medical Products, LLC or through its authorized channels of distribution;
(b) The product is altered;
(c) Any parts not supplied by Aspen Medical Products, LLC are inserted into the product; or
(d) The product is not used in accordance with the Aspen Medical Products, LLC Instructions for Use.

THE FOREGOING IS THE SOLE AND EXCLUSIVE REMEDY FOR FAILURE IN SERVICE OF, OR DEFECT IN, THE PRODUCT. ASPEN MEDICAL PRODUCTS, LLC SHALL NOT BE LIABLE UNDER THIS OR ANY IMPLIED WARRANTY FOR ANY DIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES. THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING THE WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR USE, AND ALL OBLIGATIONS OR LIABILITIES ON THE PART OF ASPEN MEDICAL PRODUCTS FOR DAMAGES ARISING OUT OF OR IN CONNECTION WITH THE USE OF THE PRODUCTS, WHICH ARE HEREBY DISCLAIMED AND EXCLUDED BY ASPEN MEDICAL PRODUCTS, LLC.

This warranty gives you specific legal rights and you may have other rights which vary from location to location.

A listing of authorized distributors in your area is available upon request.

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(d) Il prodotto non è utilizzato ai sensi delle istruzioni per l'uso fornite da Aspen Medical Products, LLC.