

ASPEN CONTOUR™ LSO COMPLETE








- LUMBAR SACRAL ORTHOSIS
- LUMBAL-SAKRALE ORTHESE
- ÓRTESIS LUMBOSACRA
- ORTOSI LOMBO-SACRALE
- ORTHÈSE LOMBO-SACRÉE
- LUMBALE SACRALE ORTHÈSE
- BEDEBNÍ KRÍŽOVÁ ORTÉZA
- LUMBAL SAKRAL ORTOSY
- ΟΣΦΥΪΚΗ ΟΡΘΩΣΗ ΟΡΤΩΣ ΟΣΤΟΥ
- LUMBAL KORSRYGGORTROE
- ORTEZA ŁĘDŹWIOWO-KRZYŻOWA
- LANNERANGAN SAKRAALINEN ORTOOSI
- LUMBAL SAKRAL ORTOS
- ОРТЕЗА ДЛЯ ПОЯСНИЧНО-КРЕСТЦОВОГО ОТДЕЛА ПОЗВОНОЧНИКА
- ÓRTESE SACRAL LOMBAR
- 腰骶矫正器

Pat. aspemp.com/patents

Aspen Medical Products, LLC
6481 Oak Canyon, Irvine, CA 92618-5202 USA
aspemp.com • +1-800-295-2776 • +1-949-681-0000

©2021 Aspen Medical Products, LLC. All Rights Reserved.
CL1306J revised 04/2021

ENGLISH

1. ASSEMBLY
For more detailed directions about the QuikDraw RAP, please visit <https://www.aspemp.com/quikdrawm-rap.html>.
A. Slide the large plastic tongue located on the posterior of the Contour into the pocket on the inside of the QuikDraw.
B. Slide the ends of the QuikDraw through the slots at the ends of the Contour. This will prevent the QuikDraw from sliding out.
C. For extra extension support, you may raise the vertical height adjustment panel in the back. Pull out the yellow rivet to disengage. When adjusted, push and snap the yellow rivet back into place.

2. APPLICACION
A. Hold the brace fully stretched behind you and center on the lower back.
B. Overlap the tapered ends of the brace so the hook grips the opposite side as tightly as possible.
C. Ends may need to be angled to accommodate different anatomies, as shown.
NOTE: You can place the tapered ends up or down to fit the anatomy of each patient.

3. TIGHTENING
A. Slip each thumb into the pull tabs. Extend the tabs outward and away from the brace.
B. When the brace is tight, press the pull tabs onto the belt to secure.

4. RAP ADJUSTMENT
If you are using a RAP, slide it as needed along the floating strap to be centered with the midline.

SIZING GUIDE		
Size	Circumference	
X-SMALL	21 – 27 in	53 – 69 cm
SMALL	26 – 32 in	66 – 81 cm
MEDIUM	31 – 37 in	79 – 94 cm
LARGE	36 – 42 in	91 – 107 cm
X-LARGE	41 – 47 in	104 – 119 cm
XX-LARGE	46 – 57 in	117 – 145 cm

5. MODIFICATION OPTIONS
This product may be substantially modified by an individual with expertise. For more information, please visit www.aspemp.com/education/product-modifications.

INTENDED USE
To provide motion restriction and stability of the lumbar-sacral spine.

INDICATIONS
Postoperative support, Deconditioned trunk musculature, Post laminectomy syndrome, Surgical Decompression procedures, Non-surgical Decompression procedures, Degenerative disc disease, Bulging or herniated disc, Non-specific low back pain, Chronic low back pain, Fracture management, Spondylolysis, Spondylolisthesis, Spinal stenosis, Ligament sprains / Muscle strains, Spondylolysis (Osteoarthritis/Degenerative Joint Disease), Facet syndrome, S1 dysfunction, Sacroiliitis, Radiculopathy, Injections, Radiofrequency Ablation, Lumbar Multifidus Muscle Atrophy.

CONTRAINDICATIONS
Hypersensitivity or allergies to any of the materials from which the brace is made. Swelling of the lymphatic tissue caused by circulatory disorders. Patients who are not allowed to wear a brace according to medical instructions.

PRECAUTIONS
READ INSTRUCTIONS BEFORE USE. Proper training in the use of this device should take place before it is applied. These directions are guidelines only and are not offered as medical recommendations. If you suffer from a serious medical condition, we strongly suggest that you consult with a licensed health care professional before using this product. Proper fitting is required for this product to be effective. Under some circumstances, this product may be prescribed by a physician. Please see the limited warranty for further information.

CARE INSTRUCTIONS
Hand wash only at 30° C. Wash with mild detergent; Air dry. Do not use bleach or other chemicals.

MATERIALS
Nylon, Cotton, Polyester, Polyurethane, HDPE, UHMWPE, Aluminum.
NOT MADE WITH NATURAL RUBBER LATEX.

DISPOSAL
Product may be safely disposed in accordance with local laws.

REPORTING
Please report any serious incident involving the use of this device to both the Competent Authority in your state and to the manufacturer (listed on this document).

LIMITED WARRANTY
Aspen Medical Products, LLC, Irvine, CA 92618, warrants to the user who originally purchases this product that it is free from defects in material and workmanship. The sole obligation of Aspen Medical Products, LLC in the event of breach of warranty shall be to repair or replace the defective product or part(s).

Aspen Medical Products, LLC shall have no obligation under this limited warranty in the event:
(a) The product was not purchased from Aspen Medical Products, LLC or through its authorized channels of distribution;
(b) The product is altered;
(c) Any parts not supplied by Aspen Medical Products, LLC are inserted into the product;
(d) The product is not used in accordance with the Aspen Medical Products, LLC Instructions for Use.

THE FOREGOING IS THE SOLE AND EXCLUSIVE REMEDY FOR FAILURE IN SERVICE OR DEFECT IN THE PRODUCT. ASPEN MEDICAL PRODUCTS, LLC SHALL NOT BE LIABLE UNDER THIS OR ANY IMPLIED WARRANTY FOR ANY DIRECT, INDIRECT, OR CONSEQUENTIAL DAMAGES. THIS WARRANTY IS VOID OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING THE WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR USE. AND ALL OBLIGATIONS OR LIABILITIES ON THE PART OF ASPEN MEDICAL PRODUCTS FOR DAMAGES ARISING OUT OF OR IN CONNECTION WITH THE USE OF THE PRODUCTS, WHICH ARE HEREBY DISCLAIMED AND EXCLUDED BY ASPEN MEDICAL PRODUCTS, LLC.

This warranty gives you specific legal rights and you may have other rights which vary from location to location.

A listing of authorized distributors in your area is available upon request.

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Aspen Medical Products, LLC non avrà alcun obbligo ai sensi della presente garanzia limitata nei seguenti casi:
(c) Il prodotto non è stato inserito parti non fornite da Aspen Medical Products, LLC; o
(d) Il prodotto non è utilizzato ai sensi delle istruzioni per l'uso fornite da Aspen Medical Products, LLC.