

ACL Knee Brace

Features and Benefits

- Polycentric hinge design to fit movement of the knee
- Prevent hyperextension via adjustable hinge ROM
- 4-leverage design for optimal stabilization
- Extension limitation at: 0°, 10°, 20°, 30°, 40°
Flexion limitation at: 45°, 60°, 75°, 90°
Immobilization limitation at: 0°, 10°, 20°, 30°
- Hypoallergenic silicone-grip strap pads for anti-migration
- Numbered straps provide desirable counter force and adjust to optimal size

Indications

Provides stability and support following injury or reconstruction of the ACL, PCL, MCL, LCL and menisci. Also used prophylactically to prevent hyperextension injuries.

CODE L1845 AND L1852 APPROVED



| Right Order # | Left Order # | Size | Thigh Circumference |
|---------------|--------------|---------|---------------------|
| 843-R | 853-L | Small | 15½ - 18½ |
| 845-R | 855-L | Medium | 18½ - 21 |
| 847-R | 857-L | Large | 21 - 23½ |
| 848-R | 858-L | X Large | 23½ - 26½ |

Note: Thigh circumference should be measured approximately 6" above the knee.