

APEX LSO/TLSO

Features

- Easy to use universal sizing system
- Adjustable pull tabs to allow optimal compression
- Breathable fabric allows air flow to keep patients cool and dry over long periods of time
- Step down 14" back panel that can be removed as patients' condition improves

Indications

Spinal Stenosis, Spondylolisthesis, Spondylolysis, Chronic Back Pain, Spinal Instability, Facet Syndrome, Degenerative Disc, Herniated Disc, Bulging Disc, Post-Operative Care

Size

Universal, Fits Waist Size 25" - 68"



Order #	Code Approved	APEX Universal
5027	L0627/L0642	APEX 627
5031	L0631/L0648	APEX 631
5037	L0637/L0650	APEX 637
5039	L0639/L0651	APEX 639
5056	L0456/L0457	APEX TLSO 456