

# Aspen<sup>®</sup> Lumbar Support

---

The Aspen Lumbar Support can help ease low back pain caused by acute flare-ups, spasms, or fatigued muscles. It also provides stability and comfort during post-operative recovery.

This inelastic brace is a low cost, medical grade treatment option that features Aspen's patented 4:1 tightening system for enhanced support. Designed to provide trunk stability and targeted compression, this low profile brace is ideal for active patients or lifestyles.

- Essential Trunk Stability
- Enhanced Level of Compression
- Tailored, Comfortable Fit
- Post Operative Recovery




# Engineered for Essential Trunk Stability and Post-Operative Recovery



## Simple Design

Comfortable for extended wear and can easily be worn under clothing.



## Independent Compression

Aspen's patented 4:1 tightening mechanism allows the patient to independently adjust the top and bottom of the brace to provide localized compression where needed.



## Integrated Anterior Panels

Anterior panels overlap to provide enhanced midline support.

## Reliable Structural Integrity

Integrated vertical stays provide structural integrity while allowing for flexibility where needed.

## Easy to Locate Pull Tabs

These ergonomic pull tabs are designed to lay flat but are easy to find by feel. Bonded and stitched materials ensure their strength and durability.



<b>X-Small</b> - 21-27 in (53-69 cm)	<b>580101</b>
<b>Small</b> - 26-31 in (66-79 cm)	<b>580102</b>
<b>Medium</b> - 31-36 in (79-91 cm)	<b>580103</b>
<b>Large</b> - 36-42 in (91-107 cm)	<b>580104</b>
<b>X-Large</b> - 41-48 in (104-122 cm)	<b>580105</b>
<b>XX-Large</b> - 43-52 in (109-132 cm)	<b>580106</b>