



COOLFORCE HIP BRACE

Features & Benefits

- Quick fitting, universal, hip range of motion orthosis
- Ambidextrous, fits left or right leg/hip
- The brace includes range of motion hinges to allow for flexion and extension control
- A quick inflation system that is easy to use, which can improve outcomes by managing pain and swelling
- Included ice pack to reduce swelling and increase recovery time

Indications

Post-Op ROM hip abduction. Indicated for use with post-op patients who have undergone a surgical procedure and need an orthosis until the component is integrated and soft tissue heals.



Order #

5525 (w/Hinges)

Size

Universal

Sugg. HCPC

L1686

THE **COOLFORCE** THERAPY LINE

This QR code can be scanned by your smart phone and directs you to the Product Detail Page.

