

CoolForce

Features and Benefits

- Liner and hand pump to force air into the brace and increase compression around the knee
- Removable ice pack designed to speed up recovery times and reduce pain
- Medial and lateral hinge stays provide excellent support
- 4-way stretch nylon fabric Piping along thigh and calf areas of support to stop flexing, Tapered in back for comfortable fit
- Hyperextension stops at 170°
- Removable hinges for custom fitting

Indications

Indicated for mild to moderate ligament joint pain, sprain, strains of the knee ligaments. Provides support, stabilization and compression of the knee joint following ligament or meniscus injuries.



Order #	Size	Code Approved
5530	Universal	N/A
5535	Universal	L1847/L1848
5545	Universal	L1810/L1812
33	Add'l Gel Pack	N/A