

# Cross-Fit Universal Hinged Knee Brace™

## Features and Benefits

- Medial and lateral hinge stays provide excellent support
- Enclosed tubular donut buttress
- Neoprene coated with 4-way stretch nylon fabric
- Piping along thigh and calf areas of support to stop flexing
- Tapered in back for comfortable fit
- Hyperextension stops at 170°
- Easy removable hinges for custom fitting
- Retail packaging available
- Universal sizing system

## Indications

Indicated for mild to moderate ligament joint pain, sprain, strains of the knee ligaments. Provides support, stabilization and compression of the knee joint following ligament or meniscus injuries.



Retail Packaging Available



Order #	Code Approved	Product
6100	L1832/L1833	Cross-Fit Universal Hinged Knee Brace