

Deluxe Shoulder Abduction System

Features and Benefits

- Lightweight, universally-sized
- Removable wedge allows shoulder positioning at 0, 15, 30, 60 and 90 degrees
- Strapping systems allows abductor to be worn on the left or right arm with minimal adjustment

Indications

Post-operative rotator cuff, elbow and shoulder injuries
Rotator cuff tear and tendinitis at the shoulder
Adhesive capsulitis at the shoulder (frozen shoulder)
Cumulative trauma/repetitive strain injury
Elbow fractures including radial head
Stroke and CVA recovery and residual hemiparesis



CODE L3960 APPROVED

Order #

Size

5100-E

Universal Fit