

# Enforcer ACL Knee Brace

## Features and Benefits

- Light weight, high strength aluminum alloy frames
- Designed to offer excellent support and protection for ACL, MCL, and other knee ligament repairs and/or instabilities
- State-of-the-art low profile ICR hinge closely mimics natural knee motion
- Customizable non-elastic straps and comfort fit condylar pads for medial-lateral stability and a superior & improved fit
- Adjustable superior and improved straps and a range of pads allow for fit customization
- Molded silicone pads offer non-slip suspension and avoid tibial pressure
- Variety of extension stops to limit range of motion
- Aluminum frames allows for additional customization by medical professionals



## Indications

Provides stability and support following injury or reconstruction of the ACL, PCL, MCL, LCL, and menisci. Can be used prophylactically to prevent hyperextension injuries.

## CODE L1845/L1852 APPROVED

Right Order #	Left Order #	Size	Knee Joint Width	Thigh Circumference
941-R	951-L	X-Small	3 - 3.5"	13.5 - 15.5"
943-R	953-L	Small	3.5 - 4"	15.5 - 18.5"
945-R	955-L	Medium	4 - 4.5"	18.5 - 21"
947-R	957-L	Large	4.5 - 5"	21 - 23.5"
948-R	958-L	X Large	5 - 5.5"	23.5 - 26.5"
949-R	959-L	2X Large	5.5 - 6"	26.5 - 29"

Note: Thigh circumference should be measured approximately 6" above the knee.