

# Ci Knee

## A Dual Approach to Ligament Instability

Our Double Upright Ligament Knee Brace is designed to address multiple types of instability to help improve confidence with increased mobility. The dual strut design helps to maintain knee stability while flexion/extension stops can be easily added to the hinge to control the allowed range of motion.

- Extension stops range from 0 - 25 degrees
- Flexion stops range from 0 (180) - 90 degrees



**CODE L1845 AND L1852 APPROVED**

SIZE	LEFT	RIGHT	CALF CIRCUMFERENCE	THIGH CIRCUMFERENCE
S	652001	652006	12-14 in (30.5-35.5 cm)	15.5-18.5 in (39-47 cm)
M	652002	652007	14-16 in (35.5-40.5 cm)	18.5-21 in (47-53.25 cm)
L	652003	652008	16-18 in (40.5-47 cm)	21-23.5 in (53.25-59.5 cm)
XL	652004	652009	18-20 in (47-50.75 cm)	23.5-26.5 in (59.5-67.25 cm)

# A dual approach to ligament instability.



## Easy Application and Removal

Proprietary one-step closure system for easy donning and doffing.



## Dual Strut Design

Provides increased stability.



## Advanced Aerospace Grade Aluminum

Provides stability and support for ligament injuries such as ACL, PCL and combined instabilities.



## Flexion/Extension Adjustments

Easily control the allowed range of motion.



Ci Knee, Strap Accessory Kit – S	652120	Ci Knee, Strap Accessory Kit – M	652121
Ci Knee, Strap Accessory Kit – L	652122	Ci Knee, Strap Accessory Kit – XL	652123
Ci Knee, Liner Accessory Kit – S	652125	Ci Knee, Liner Accessory Kit – M	652126
Ci Knee, Liner Accessory Kit – L	652127	Ci Knee, Liner Accessory Kit – XL	652128
Ci Knee, Extension/Flexion Stop Frame	652130	Aspen Quick Release Buckle Kit	651130