

Horizon™ PRO Lower Spine Braces

Backed by clinical research and 10 years of trusted technology, the Horizon PRO provides an updated approach to pain management. Its inelastic design combined with targeted compression offers immediate relief from symptoms caused by injury, muscle fatigue, and spasms.

One Size Adjustable

Comfortably fits waists ranging from 24 to 50 inches and up to 70 inches with the extension panels. Simple pull-through mechanism offers multiple sizes in one brace, enabling each brace to comfortably fit various patient shapes and sizes. One size adjustable means the correct size is always on hand, minimizing inventory requirements and significantly reducing costs.



NEW!



Horizon PRO 464



Horizon PRO 456

Manage symptoms of low back pain and add trunk stability with market-leading technology.

Integrated Anterior Panels

Anterior panels overlap to provide enhanced midline support. Featuring Aspen's innovative FlexTabs™, these panels disperse pressure and naturally conform to the anatomy for added comfort.

Low Profile Design with Reliable Structural Integrity

The Horizon PRO 639 features a tall, rigid anterior panel and overlapping lateral/anterior panels that provide flexion and lateral motion restriction without compromising comfort. The xiphoid extension can be angled as needed or removed for shorter patients.




Ergonomic Pull Tabs

Comfortable and easy to use for patients to achieve their desired level of compression.



Easy Access to Panels

Improve customization in the clinic.



Custom Support

When tightening, the firm but flexible back panel naturally conforms to the patient's lordosis to provide added structural support.

Integrated Pockets

Provide easier application and ensure a snug fit.




Height Adjusts for Every Body

Telescoping shoulder straps on the Horizon PRO 456 provide 13 inches of height adjustment to fit varying patient anatomies. In addition, circumferential adjustability enables a customized fit for every body.

Independent Tightening Mechanism

Aspen's patented 4:1 mechanically advantaged SlickTrack™ System features an independent upper and lower tightening mechanism, providing direct compression where needed to support tight muscles.



The Right Brace for Multiple Needs

Regardless of which configuration a patient requires initially, all Horizon PRO braces can be stepped down to lowest-profile option for the treatment of reoccurring ailments.

← FITS WAIST CIRCUMFERENCE RANGING FROM 24-50 IN (61-127 CM) AND UP TO 70 IN (178 CM) WITH HORIZON EXTENSION PANEL →

Horizon™ PRO 464 TLSO - Code L0464 Approved	994760	Horizon™ PRO 631 LSO - Code L0631 & L0648 Approved	994720
Horizon™ PRO 456 TLSO - Code L0456 & L0457 Approved	994740	Horizon™ PRO 627 Lumbar - Code L0627 & L0642 Approved	994710
Horizon™ PRO 639 LSO - Code L0639 & L0651 Approved	994735	Horizon™ Extension Panels - 20 in (51 cm)	993745
Horizon™ PRO 637 LSO - Code L0637 & L0650 Approved	994730	Therapy Pack	990005