

# FixStep Walker

**Fixed Strut Walker available in both standard and short sizes with either high volume pump air liners or standard liners**

- Excellent rocker sole profile for improved gait
- Two-part injection molded sole unit
- Improved shock absorbing - less impact during walking
- Flared & deeper heel section to retain foot and reduce any pressure on Achilles area
- High quality hook/loop straps - easy for patients to open
- Strap lock system avoids straps rotating
- Molded high volumn pump on air liner for faster inflation and easier patient use
- Comfortable detachable liners with replacements available
- Reinforced molded plastic uprights - giving the FixStep a metal free construction and making the boots x-ray safe



**Air Liner Options**



**Easy Open Straps**



**Rocker Sole**

							
Size	FixStep with Standard Liner L4386/L4387 Approved	Replacement Standard Liner for FixStep	Short FixStep with Standard Liner L4386/L4387 Approved	Replacement Standard Liner for Short FixStep	FixStep with Air Liner L4360/L4361 Approved	Short FixStep with Air Liner L4360/L4361 Approved	Shoe Size
XS	FXP/02	FXPL/02	FXPS/02	FXPSL/02	FXPA/02	FXPSA/02	2.5 - 4
S	FXP/03	FXPL/03	FXPS/03	FXPSL/03	FXPA/03	FXPSA/03	5 - 6.5
M	FXP/05	FXPL/05	FXPS/05	FXPSL/05	FXPA/05	FXPSA/05	7 - 10
L	FXP/07	FXPL/07	FXPS/07	FXPSL/07	FXPA/07	FXPSA/07	10.5 - 12
XL	FXP/08	FXPL/08	FXPS/08	FXPSL/08	FXPA/08	FXPSA/08	12.5 - 14.5