

DuAlign OA

A Dual Approach to OA

The DuAlign OA Knee Brace is designed to help increase mobility, comfort, and confidence by addressing pain from mild to severe OA. The dual strut design helps maintain knee stability while active. Dual sided varus/valgus adjustment allows for more precise control of offloading and positioning. Extension stops can be easily added to the hinge to change the allowable range of motion.

- Extension stops range from 0 - 25 degrees
- Flexion stops range from 0 (180) - 90 degrees



CODE L1845 AND L1852 APPROVED

SIZE	LEFT	RIGHT	CALF CIRCUMFERENCE	THIGH CIRCUMFERENCE
S	652011	652016	12-14 in (30.5-35.5 cm)	15.5-18.5 in (39-47 cm)
M	652012	652017	14-16 in (35.5-40.5 cm)	18.5-21 in (47-53.25 cm)
L	652013	652018	16-18 in (40.5-47 cm)	21-23.5 in (53.25-59.5 cm)
XL	652014	652019	18-20 in (47-50.75 cm)	23.5-26.5 in (59.5-67.25 cm)

A dual approach to OA.



Easy Application and Removal

Proprietary one-step closure system for easy donning and doffing.



Flexion/Extension Adjustments

Easily control the allowable range of motion.



Two Point Hinge Adjustment

2 Hinge adjustment points, one on either strut, allow for even pressure distribution, improved comfort, and more precise medial or lateral correction.



Dual Strut Design

Provides increased stability and security.



Strap Accessory Kit – S	652120	Strap Accessory Kit – M	652121
Strap Accessory Kit – L	652122	Strap Accessory Kit – XL	652123
Liner Accessory Kit – S	652125	Liner Accessory Kit – M	652126
Liner Accessory Kit – L	652127	Liner Accessory Kit – XL	652128
Extension/Flexion Stop Frame	652130	Aspen Quick Release Buckle Kit	651130