

# Universal Wrist Brace with Thumb Spica

## Features and Benefits

- Universal sizing
- Adjustable dorsal stay
- Single strap design
- Constructed of lightweight, breathable material
- QR Code tag for patient education
- Removable and malleable aluminum palmer stay for correct fitting

## Indications

Wrist immobilization, De Quervain's Syndrome, wrist and thumb sprains and strains

**CODE L3807/L3809 APPROVED**



Retail Packaging  
Available



**Order #**

**Size**

**170-R**

Right Universal

**180-L**

Left Universal