





## Giving Experts Modification Options

# MOLD BEND TRIM


HIGH QUALITY MATERIALS



Easy Access



Mold Panels



Trim Sure Slot™



Trim Panels



Bend TLSO T-Bar



Bend Support Bar



### DESIGNED FOR MODIFICATION

Aspen® Spinal Braces provide unique opportunities to trim, bend, mold, assemble and provide individualized customization. The various structural elements that comprise each brace are accessible to the provider. This allows the practitioner the ability to remove, heat mold, bend, and trim these elements in order to modify the brace for individual patient needs.

### SPECIFIC PATIENT INDICATIONS

Aspen offers numerous options to treat specific patient indications and aid with treatment of different diagnoses. Each practitioner

is able to modify the brace to adapt to different bracing treatment philosophies or physician direction. Depending on the clinical justification, Aspen Braces can be modified to influence the patient's spinal alignment and accommodate for additional medical interventions. Aspen recognizes the importance of the clinical care required to meet the medical needs of each patient. With this in mind, Aspen's goal is to develop braces that can be significantly modified as needed, leading to improved outcomes and increased patient compliance.

### ASPEN QUALITY OUTPERFORMS

Not all braces perform the same. Aspen's

Spinal Braces are made of the highest quality materials such as nylon, high-density and low-density polyethylene plastics and non-ferrous metals. These durable and resilient materials not only provide the best support, they allow for significant modification. Heating, bending and molding lower quality materials can cause cracking, buckling, and breakage. Aspen offers quality products you can count on.

### CLINICAL PROVIDER DOCUMENTATION

Aspen is committed to providing resources for practitioners regarding modification documentation. For more information, visit us at [aspenmp.com](http://aspenmp.com).



Life Changing Spinal Innovation



The Aspen Advantage Symbol

An innovative feature available only from Aspen