

# Aspen Hinged Knee

Maximum Knee Support In A Low-Profile Design

XS	648001	L	648004
S	648002	XL	648005
M	648003	XXL	648006

The Aspen Hinged Knee sized brace features a premium AirX material that provides increased compression and support with a high degree of breathability and moisture wicking. The comfort provided by the AirX material also allows for extended wear and improves patient compliance.

**CODE L1820 APPROVED**




#### Aspen Hinged Knee Universal

Provides adjustability to fit most patients. Made from a lightweight neoprene material, this option allows for optimal compression and comfort.

**Aspen Hinged Knee Universal**  
Code L1820 Approved

**648011**





#### Comfort

Advanced AirX material provides ideal compression with a high degree of breathability and moisture wicking.



#### Hinge

Polycentric hinge for peak performance.



#### Pull Tabs

Molded pull tabs for easy adjustment.



#### Condyle Pads

Integrated condyle pads for added comfort.



## Aspen Knee Suspension Wrap

The Aspen Knee Suspension Wrap can be used in conjunction with the Aspen Hinged Knee for patients needing additional support.

#### Accessory

**Aspen Knee Suspension Wrap**  
Code L2397 Approved

**651111**

<b>Aspen Hinged Knee XS</b> - Code L1820 Approved - 11-15 in (27.94-38.1 cm)	<b>648001</b>
<b>Aspen Hinged Knee S</b> - Code L1820 Approved - 15-18 in (38.1-45.72 cm)	<b>648002</b>
<b>Aspen Hinged Knee M</b> - Code L1820 Approved - 18-21 in (45.72-53.34 cm)	<b>648003</b>
<b>Aspen Hinged Knee L</b> - Code L1820 Approved - 21-24 in (53.34-60.96 cm)	<b>648004</b>
<b>Aspen Hinged Knee XL</b> - Code L1820 Approved - 24-27 in (60.96-68.58 cm)	<b>648005</b>
<b>Aspen Hinged Knee XXL</b> - Code L1820 Approved - 27-35 in (68.58-88.9 cm)	<b>648006</b>
<b>Aspen Hinged Knee Universal</b> - Code L1820 Approved	<b>648011</b>