




## Aspen Quality Sacral Stabilization

The new Evergreen™ SI delivers the quality of design and durability that you've come to expect from Aspen.

Code L0621 Approved



## Low Profile and Lightweight


Durable two-ply construction with breathable mesh fabric makes the Evergreen SI comfortable for extended wear.

## Bi-Lateral Pull

Two easy-to-locate pull tabs provide maximum compression without the spinal rotation that can occur with a single pull tab system.

## Effective Compression

Utilizes Aspen's proven SlickTrack™ mechanically advantaged tightening system. Easy to pull regardless of patient strength.



## Non-Slip Strips

Silicone strips on the inside of the brace hold fast and prevent migration especially while moving and sitting.



## Product Information

Evergreen SI, Small	36-2102
Evergreen SI, Medium	36-2104
Evergreen SI, Large	36-2106