



Raising the Level of Patient Care

Six Sizes in One

The Vista Cervical Collar from Aspen improves patient care while saving time, money and storage space. Built on the design of the trusted Aspen Collar, the Vista Collar provides safe and effective motion restriction.¹ With its innovative height adjustment technology, the Vista Collar is really six sizes in one. The right size is always at hand, reducing storage and inventory costs, while improving patient care.

1. Comparative Motion Allowed, Vista vs. Aspen, Aspen Medical Products in-house testing, December 2006

Code L0172 Approved


Eliminate Sizing

With six sizes in one, the Vista Collar eliminates the need for pre-sizing methods of any kind.



Unparalleled Patient Access

The extra large tracheal aperture provides functional access for airway management and other essential procedures.



Dial Height Adjustment

Simply pull and turn the Height Adjustment Dial to select one of six Vista height settings.



Product Information

Vista Collar ————— **984000**
Vista Collar Set ————— **984002**

Inventory Management

By incorporating six sizes in one collar, the Vista Collar eliminates the cost of storing multiple sizes, frees up valuable inventory space, and eliminates SKUs. The Vista Collar reduces the waste associated with sizing errors and eliminates the need to change collars as patients become ambulatory.



Life Changing Spinal Innovation



The Aspen Advantage Symbol

An innovative feature available only from Aspen

vistaicu

BACK PANEL

Inspired by ICU Nurses

Maximum Skin Care and Comfort

Long known and respected for cervical care products that promote quality skin care without sacrificing comfort and motion restriction, Aspen is proud to introduce our newest innovation in spine care. The Vista ICU Back Panel is designed specifically for patients at risk of occipital breakdown due to extended periods lying supine.



Pressure Reducing Pad

The Vista ICU Back Panel is designed to disperse pressure, reducing the possibility of skin breakdown. Proven Aspen Pads^{1,2} complete the system by wicking moisture away from the skin.

1. Powers, J. *A Multidisciplinary Approach to Occipital Pressure Ulcers Related to Cervical Collars.* *Journal of Nursing Quality*, 1997, Vol. 2, No. 1
2. *Absorption Comparison Test, Orange County Materials Test Laboratories, Anaheim CA, May 28, 2002*

*Vista ICU Back Panel and Vista Collar sold separately.



Extended Back Panel Height


Prevents bed shear when rotating or repositioning patients.

Unique Dual-Density Foam

Low density white ether foam cradles patient's head for maximum pressure distribution. Higher density grey polyurethane foam provides added support and structure.

Open by Design

Open design promotes airflow for ventilation without sacrificing support.



Effective Motion Restriction

The Vista ICU Back Panel utilizes the standard Vista Back Panel providing the same superior motion restriction. Back panel extension can be removed once patient becomes ambulatory.

Product Information

Vista ICU Back Panel ————— **984600**
 Vista ICU Back Panel Replacement Pads ————— **984620**